


Biotin Rabbit Anti-Human Cd40 Monoclonal Antibody For Fc - P25942

Item Number: CM0008439



Introduction

Biotin-conjugated rabbit monoclonal antibody against human CD40, validated for flow cytometry. Reacts with human CD40, a TNF receptor superfamily member involved in immune response. Suitable for diagnostic and research use.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Human CD40 mAb
Remarks/Alias	p50; Bp50; CDW40; TNFRSF5
Species	Human
Gene ID	958
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-193 of human CD40 (NP_001241.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P25942
Protein Weight	22kDa/31kDa
Shipping	Ice bag