

Biotin Rabbit Anti-Human Cd45Ra Monoclonal Antibody For Flow Cytometry - P08575-8

Item Number: CM0010586

Introduction

Biotinylated rabbit monoclonal antibody targeting human CD45RA (Isoform 6 of PTPRC). Validated for flow cytometry (FC) applications. Ideal for leukocyte subtyping in research and IVD assay development. Molecular weight ~138 kDa. Reacts with human.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Human CD45RA mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
GeneID (Human)	5788
GeneID	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-482 of human CD45RA (P08575-8)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08575-8
Protein Weight	138kDa
Shipping	Ice bag