

Biotin Rabbit Anti-Human/Monkey Cd14 Monoclonal Antibody For Fc - P08571

Item Number: CM0005390



Introduction

Biotin-labeled rabbit monoclonal antibody targeting human and cynomolgus CD14, validated for flow cytometry. Recombinant immunogen; ensures reactivity with monocytes and macrophages; ideal for innate immunity and LPS signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Human/Monkey CD14 mAb
Remarks/Alias	CD14
Species	Human
GeneID (Human)	929
GeneID	929
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-352 of human CD14 (NP_000582.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08571
Protein Weight	40 kDa
Shipping	Ice bag