

Biotin Rabbit Anti-Human/Monkey Cd16 Monoclonal Antibody For Fc - P08637 / O75015

Item Number: CM0011673



Introduction

Biotinylated rabbit monoclonal antibody against human/monkey CD16 (FCGR3A/FCGR3B) for flow cytometry. Cross-reacts with human and cynomolgus monkey, ideal for NK cell and ADCC research.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Human/Monkey CD16 mAb
Remarks/Alias	CD16; FCG3; CD16A; FCGR3; IGFR3; IMD20; FCR-10; FCRIII; CD16-II; FCGRIII; FCRIIIA; FcGRIIIA
Species	Human
Gene ID (Human)	2214/2215
Gene ID	2214/2215
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-208 of human FCGR3A/FCGR3B (P08637/O75015)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08637/O75015
Protein Weight	29KDa/26KDa
Shipping	Ice bag