

Biotin Rabbit Anti-Human/Monkey Cd3 Monoclonal Antibody For Flow Cytometry - P07766

Item Number: CM0010608



Introduction

Biotin-conjugated rabbit monoclonal antibody targeting human and cynomolgus monkey CD3 epsilon. Validated for flow cytometry, this reagent supports T-cell profiling in immuno-oncology and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Human/Monkey CD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
Gene ID (Human)	916
Gene ID	916
Immunogen	Recombinant protein (amino acids 1-126 of human CD3, NP_000724.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag