

Biotin Rabbit Anti-Human/Monkey Cd3 Monoclonal Antibody For Flow Cytometry - P07766

Item Number: CM0008181



Introduction

Biotin-conjugated rabbit monoclonal antibody against human and cynomolgus CD3 epsilon chain, validated for flow cytometry. Targets the TCR-CD3 complex, essential for T-cell activation and development studies.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Human/Monkey CD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
Gene ID (Human)	916
Gene ID	916
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag