

Biotin Rabbit Anti-Human/Monkey Hla-Dr Mab For Fc - P01903

Item Number: CM0012367



Introduction

Biotin-conjugated rabbit monoclonal antibody against HLA-DR alpha chain, validated for flow cytometry on human and cynomolgus samples. Ideal for immune cell phenotyping and MHC class II antigen presentation studies.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Human/Monkey HLA-DR mAb
Remarks/Alias	HLA-DRA1
Species	Human
Gene ID (Human)	3122
Gene ID	3122
Immunogen	Recombinant protein of human HLA-DRA
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01903
Protein Weight	29kDa
Shipping	Ice bag