

Biotin Rabbit Anti-Mouse Cd19 Monoclonal Antibody For Flow Cytometry - P25918

Item Number: CM0009868

Introduction

Biotin-conjugated rabbit monoclonal antibody specific for mouse CD19 (UniProt P25918), validated for flow cytometry. Enables detection of B lymphocytes in immunophenotyping and B-cell development studies.

[Learn More](#)

Parameter	Description
Product Name	Biotin Rabbit anti-Mouse CD19 mAb
Species	Mouse
Gene ID (Human)	930
Gene ID	12478
Immunogen	Recombinant protein of mouse CD19
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P25918
Protein Weight	60kDa
Shipping	Ice bag