

Biotin Rabbit Anti-Mouse Cd25 Monoclonal Antibody For Flow Cytometry - P01590

Item Number: CM0006151

Introduction

Biotin-conjugated rabbit monoclonal antibody targeting mouse CD25 (IL2RA) for flow cytometry (FC). CD25 is the interleukin-2 receptor alpha chain, key in immune tolerance and Treg function. Suitable for mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Mouse CD25 mAb
Remarks/Alias	CD25; IL2r; Ly-43
Species	Mouse
Gene ID (Human)	3559
Gene ID	16184
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-236 of mouse CD25 (NP_032393.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01590
Protein Weight	30kDa
Shipping	Ice bag