

Biotin Rabbit Anti-Mouse Cd28 Monoclonal Antibody For Flow Cytometry - P31041

Item Number: CM0012572



Introduction

Biotin-conjugated rabbit monoclonal antibody against mouse CD28, validated for flow cytometry. Ideal for T-cell activation and costimulation studies in immunological research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Mouse CD28 mAb
Remarks/Alias	Cd28
Species	Mouse
GeneID (Human)	940
GeneID (Mouse)	12487
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-149 of mouse CD28 (NP_031668.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P31041
Protein Weight	25kDa
Shipping	Ice bag