

Biotin Rabbit Anti-Mouse Cd4 Mab - P06332

Item Number: CM0011372

Introduction

Biotin-conjugated rabbit monoclonal antibody targeting mouse CD4 (L3T4, Ly-4) for flow cytometry. Immunogen covers extracellular domain; validated for mouse. Ideal for T-helper cell and immune response studies.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Mouse CD4 mAb
Remarks/Alias	L3T4; Ly-4
Species	Mouse
Gene ID (Human)	920
Gene ID	12504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-394 of mouse CD4 (NP_038516.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06332
Protein Weight	51 kDa
Shipping	Ice bag