

Biotin Rabbit Anti-Mouse Cd4 Monoclonal Antibody For Flow Cytometry - P06332

Item Number: CM0005197

Introduction

Biotinylated rabbit monoclonal antibody targeting mouse CD4 (L3T4, Ly-4). Validated for flow cytometry on mouse samples. Coreceptor for MHC class II essential for T-helper cell function. Swiss-Prot: P06332.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Mouse CD4 mAb
Remarks/Alias	L3T4; Ly-4
Species	Mouse
GeneID (human)	920
GeneID	12504
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 27-394 of mouse CD4 (NP_038516.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06332
Protein Weight	51kDa
Shipping	Ice bag