

Biotin Rabbit Anti-Mouse Cd45 Monoclonal Antibody For Flow Cytometry - P06800

Item Number: CM0009244

Introduction

Biotin-conjugated rabbit monoclonal antibody against mouse CD45 (PTPRC, L-CA, B220) for flow cytometry. Suitable for leukocyte common antigen detection in mouse samples, validated for FC.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Mouse CD45 mAb
Remarks/Alias	loc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyt-4
Species	Mouse
GenelD (Human)	5788
GenelD	19264
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-566 of mouse CD45 (NP_001104786.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800
Protein Weight	130 kDa/140 kDa/144 kDa
Shipping	Store at -20°C