

Biotin Rabbit Anti-Mouse Cd45 Monoclonal Antibody For Flow Cytometry - P06800

Item Number: CM0008686

Introduction

Biotinylated rabbit monoclonal antibody against mouse CD45 (P06800), validated for flow cytometry. Recombinant immunogen (aa 26-566), detects mouse CD45 with multiple isoforms (130-144 kDa) for immunology research.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Mouse CD45 mAb
Remarks/Alias	loc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyt-4
Species	Mouse
GeneID (Human)	5788
GeneID	19264
Immunogen	Recombinant fusion protein containing amino acids 26-566 of mouse CD45 (NP_001104786.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800
Protein Weight	130kDa/140kDa/144kDa
Shipping	Ice bag