

Biotin Rabbit Anti-Mouse Cd5 Monoclonal Antibody For Flow Cytometry - P13379

Item Number: CM0011224

Introduction

Biotin-conjugated rabbit monoclonal antibody against mouse CD5 (Ly-1). Designed for flow cytometry (FC) detection of T-cell surface glycoprotein in mouse samples. Validated for mouse cross-reactivity. 54 kDa target.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Mouse CD5 mAb
Remarks/Alias	Ly-1; Ly-A; Ly-12; Lyt-1
Species	Mouse
Gene ID (human)	921
Gene ID	12507
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-371 of mouse CD5 (NP_031676.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P13379
Protein Weight	54kDa
Shipping	Ice bag