


Biotin Rabbit Anti-Mouse Cd8A Monoclonal Antibody For Fc - P01731

Item Number: CM0010692



Introduction

Biotin-conjugated rabbit monoclonal antibody targeting mouse CD8a (Ly-2) for flow cytometry. Recognizes CD8 alpha chain, a 27 kDa T-cell coreceptor critical for MHC class I interaction and CTL-mediated immunity.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Mouse CD8a mAb
Remarks/Alias	Ly-2; Ly-B; Ly-35; Lyt-2
Species	Mouse
Gene ID (Human)	925
Gene ID	12525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-196 of mouse CD8a (NP_001074579.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01731
Protein Weight	27kDa
Shipping	Ice bag