

Biotin Rabbit Anti-Mouse Cd8A Monoclonal Antibody For Fc - P01731

Item Number: CM0009201



Introduction

Biotin-conjugated rabbit monoclonal antibody against mouse CD8a for flow cytometry. Specific for mouse T-cell surface glycoprotein CD8 alpha chain (MW 27kDa). Ideal for immune response studies.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Mouse CD8a mAb
Remarks/Alias	Ly-2; Ly-B; Ly-35; Lyt-2
Species	Mouse
GeneID (Human)	925
GeneID	12525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-196 of mouse CD8a (NP_001074579.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01731
Protein Weight	27kDa
Shipping	Ice bag