

Biotin Rabbit Anti-Mouse IgA (Fc) Mab

Item Number: CM0008678

Introduction

Biotinylated rabbit anti-mouse IgA (Fc) monoclonal antibody for intracellular flow cytometry. Specific to mouse IgA Fc region, ensuring minimal interference. Ideal for immunology and mucosal immunity research.

[Learn More](#)

Parameter	Specification
Product Name	Biotin Rabbit anti-Mouse IgA (Fc) mAb
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 99-344 of mouse IgA.
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow cytometry (intracellular)
Cross Reactivity	Mouse
Shipping	Ice bag