

Biotin Rabbit Anti-Mouse IgM Monoclonal Antibody For Flow Cytometry - P01872

Item Number: CM0004827



Introduction

Biotin-conjugated rabbit monoclonal antibody against mouse IgM heavy chain constant region, specifically validated for flow cytometry. Detects both membrane-bound BCR and secreted IgM on mouse B cells.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Mouse IgM mAb
Remarks/Alias	Igm; muH; Igh6; Igh-6; Igh-M
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 105-454 of mouse IgM (P01872)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01872
Protein Weight	50kDa
Shipping	Ice bag