

Biotin Rabbit Anti-Mouse Integrin Alpha 2/Cd49B Monoclonal Antibody For Fc - Q62469

Item Number: CM0005661

Introduction

Biotin-conjugated rabbit monoclonal antibody targeting mouse integrin alpha-2 (CD49b/DX5/GPIa). Validated for flow cytometry. Suitable for collagen receptor and cell adhesion studies. Predicted molecular weight 129 kDa.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Mouse Integrin alpha 2/CD49b mAb
Remarks/Alias	DX5; GPIa; CD49B
Species	Mouse
GenelD (human)	3673
GenelD	16398
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 168-366 of mouse Integrin alpha 2/CD49b (NP_032422.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q62469
Protein Weight	129 kDa
Shipping	Ice bag