

Biotin Rabbit Monoclonal Antibody For Biotin Detection In Wb And Elisa Applications

Item Number: CM0000014

Introduction

This biotin rabbit monoclonal antibody targets biotin and is validated for WB and ELISA applications. It is suitable for biotin labeling research, detection of biotinylated molecules, and quality control in diagnostic development.

[Learn More](#)

Parameter	Specification
Product Name	Biotin Rabbit mAb
Species	Human
Immunogen	Recombinant protein, Biotin
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, ELISA
Cross Reactivity	Other (Wide Range Predicted)
Shipping	Ice bag