

Biotin Rabbit Polyclonal Antibody

Item Number: CM0021067

Introduction

Biotin Rabbit Polyclonal Antibody (pAb) for WB detection. Species-independent cross-reactivity enables broad use in biotinylation assays and biotin-labeled protein analysis. Supplied with ice bag shipping.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit pAb
Immunogen	Chemical compounds
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB
Cross Reactivity	Species independent
Shipping	Ice bag