

Biotin Rat Anti-Mouse Cd3 Monoclonal Antibody For Fc - P04235 / P22646

Item Number: CM0004415

Introduction

Biotin-conjugated rat monoclonal antibody against mouse CD3 delta and epsilon chains. Validated for flow cytometry, ideal for studying T-cell activation, thymocyte differentiation, and adaptive immunity.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rat anti-Mouse CD3 mAb
Remarks/Alias	T3d; CD3; T3e; CD3epsilon
Species	Mouse
GenelD (Human)	916/915
GenelD	12500/12501
Immunogen	A synthesized peptide derived from mouse CD3.
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P04235/P22646
Protein Weight	19kDa/21kDa
Shipping	Ice bag