

Biotin Rat Anti-Mouse Ter-119 Mab - Ter-119

Item Number: CM0006219

Introduction

Biotin-conjugated rat monoclonal antibody targeting mouse TER-119 antigen for flow cytometry (FC) applications. Specific to mouse samples, suitable for research use in erythroid lineage studies.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rat anti-Mouse TER-119 mAb
Remarks/Alias	Ter119; TER-119
Species	Mouse
Gene ID	104231
Immunogen	Recombinant protein of mouse TER-119
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
Shipping	Ice bag