

Bmpr1B Rabbit Pab - 000238

Item Number: CM0019402



Introduction

BMPR1B Rabbit pAb is a polyclonal antibody validated for IF-P, IHC-P, and ELISA. It detects human, mouse, and rat BMPR1B, a receptor for BMP7 and GDF5 involved in chondrocyte differentiation.

[Learn More](#)

Parameter	Value
Product Name	BMPR1B Rabbit pAb
Remarks/Alias	ALK6; AMD3; AMDD; BDA2; ALK-6; BDA1D; CDw293; BMPR1B
Species	Human
Gene ID (Human)	658
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00238
Protein Weight	57kDa
Shipping	Ice bag