

Bmpr2 Rabbit Pab For Wb, If/Icc, Elisa - Q13873

Item Number: CM0029940



Introduction

BMPR2 Rabbit pAb (Q13873) is a polyclonal antibody targeting bone morphogenetic protein receptor type-2 for WB, IF/ICC, ELISA. With cross-reactivity to human, mouse, rat, it suits pulmonary hypertension and BMP signaling studies.

[Learn More](#)

Parameter	Value
Product Name	BMPR2 Rabbit pAb
Remarks/Alias	BMR2; PPH1; BMPR3; BRK-3; POVD1; T-ALK; BMPR-II; BMPR2
Species	Human
GeneID (Human)	659
GeneID	659
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-150 of human BMPR2 (NP_001195.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13873
Protein Weight	115kDa
Shipping	Ice bag