

Bob1 (Pou2Af1) Rabbit Polyclonal Antibody For Wb, Elisa - Q16633

Item Number: CM0020803



Introduction

High specificity BOB1 (POU2AF1) rabbit polyclonal antibody validated for Western blotting (WB) and ELISA. Targets human BOB1 transcription coactivator involved in B-cell immune response, ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	BOB1 Rabbit pAb
Remarks/Alias	POU2AF1; BOB1; OBF-1; OBF1; OCAB; POU class 2 associating factor 1
Species	Human
GeneID	5450
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-250 of human BOB1 (NP_006226.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q16633
Protein Weight	27kDa
Shipping	Ice bag