

Brd3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15059

Item Number: CM0019116



Introduction

Rabbit polyclonal antibody targeting human BRD3, validated for Western blot, immunohistochemistry-paraffin, and ELISA. Detects human and rat BRD3. Ideal for chromatin and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	BRD3 Rabbit pAb
Remarks/Alias	ORFX; FSHRG2; RING3L; BRD3
Species	Human
Gene ID (Human)	8019
Gene ID	8019
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-305 of human BRD3 (NP_031397.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q15059
Protein Weight	80kDa
Shipping	Ice bag