


Brd4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O60885

Item Number: CM0021369



Introduction

Rabbit polyclonal antibody against human BRD4 (bromodomain-containing protein 4), suitable for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, and rat BRD4. Useful in epigenetic and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	BRD4 Rabbit pAb
Remarks/Alias	CAP; MCAP; HUNK1; HUNKI; FSHRG4; BRD4
Species	Human
Gene ID	23476
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 700-794 of human BRD4 (NP_001317313.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60885
Protein Weight	152 kDa
Shipping	Ice bag