


Brd7 Rabbit Pab - Q9Npi1

Item Number: CM0020214



Introduction

Anti-BRD7 rabbit polyclonal antibody for WB and ELISA; reacts with human, mouse, and rat. BRD7 is a bromodomain protein involved in chromatin remodeling, Wnt signaling, and tumor suppression. Suitable for transcription and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	BRD7 Rabbit pAb
Remarks/Alias	BP75; NAG4; CELTIX1; SMARCI1; BRD7
Species	Human
Gene ID (Human)	29117
Gene ID	29117
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human BRD7 (NP_037395.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPI1
Protein Weight	74kDa
Shipping	Ice bag