


## Bre Rabbit Polyclonal Antibody - Q9Nxr7

Item Number: CM0016366



### Introduction

BRE Rabbit Polyclonal Antibody validated for WB, IF/ICC, IF-P, ELISA. Cross-reacts with Human, Mouse, Rat. Immunogen: recombinant protein (aa 51-229). Used in DNA damage repair and apoptosis research.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | BRE Rabbit pAb                                                                                                  |
| Remarks/Alias    | BRE; BRCC4; BRCC45                                                                                              |
| Species          | Human                                                                                                           |
| GeneID (human)   | 9577                                                                                                            |
| GeneID           | 9577                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 51-229 of human BRE (NP_954661.1) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, IF/ICC, IF-P, ELISA                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                               |
| SWISS            | Q9NXR7                                                                                                          |
| Protein Weight   | 44kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |