


Bscl2 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - Q96G97

Item Number: CM0016638



Introduction

High-quality rabbit polyclonal antibody against human BSCL2 (seipin). Suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for lipid droplet biology and adipogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	BSCL2 Rabbit pAb
Remarks/Alias	HMN5; PELD; HMN5C; SPG17; GNG3LG; BSCL2
Species	Human
GeneID (Human)	26580
GeneID	26580
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 265-363 of human BSCL2 (NP_116056.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96G97
Protein Weight	44kDa
Shipping	Ice bag