


Bst1 Rabbit Polyclonal Antibody - Q10588

Item Number: CM0025367



Introduction

BST1 rabbit polyclonal antibody (CD157, cADPR2) validated for Western blot, immunofluorescence, and ELISA. Reacts with human and mouse samples. Suitable for studying cyclic ADP-ribose signaling and pre-B-cell growth.

[Learn More](#)

Parameter	Value
Product Name	BST1 Rabbit pAb
Remarks/Alias	CD157; cADPR2; BST1
Species	Human
GeneID (Human)	683
GeneID	683
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 29-293 of human BST1 (NP_004325.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q10588
Protein Weight	36kDa
Shipping	Ice bag