


Btk Rabbit Polyclonal Antibody For Wb, Elisa - Q06187

Item Number: CM0013347



Introduction

BTK rabbit polyclonal antibody validated for WB and ELISA, with cross-reactivity to mouse and rat. Target tyrosine-protein kinase BTK (Q06187), a key B-cell signaling molecule, for immunology research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	BTK Rabbit pAb
Remarks/Alias	AT; ATK; BPK; XLA; IMD1; AGMX1; IGHD3; PSCTK1; BTK
Species	Human
GeneID (human)	695
GeneID	695
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 173-272 of human BTK (NP_000052.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q06187
Protein Weight	76kDa
Shipping	Ice bag