


C-Abl (Abl1) Rabbit Monoclonal Antibody For Wb, Elisa - P00519

Item Number: CM0007670



Introduction

Recombinant rabbit monoclonal antibody against human c-Abl (ABL1) for Western blot and ELISA. Detects the 123 kDa non-receptor tyrosine kinase involved in cell growth, cytoskeleton remodeling, and apoptosis. Suitable for cancer and signaling research.

[Learn More](#)

Parameter	Value
Product Name	c-Abl Rabbit mAb
Remarks/Alias	ABL; JTK7; p150; c-ABL; v-abl; CHDSKM; c-ABL1; BCR-ABL; bcr/abl; c-Abl
Species	Human
GeneID (Human)	25
GeneID	25
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1000-1130 of human c-Abl (NP_005148.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P00519
Protein Weight	123kDa
Shipping	Ice bag