

C-Maf Rabbit Monoclonal Antibody For Wb And Elisa - 075444

Item Number: CM0008176

Introduction

Rabbit monoclonal c-Maf antibody validated for WB, ELISA with human, mouse, rat reactivity. Targets transcription factor Maf (38kDa). Ideal for lens development, Th2 cell, and oncogene/tumor suppressor research.

[Learn More](#)

Parameter	Value
Product Name	c-Maf Rabbit mAb
Remarks/Alias	CCA4; AYGRP; c-MAF; CTRCT21; c-Maf
Species	Human
GeneID (human)	4094
GeneID	4094
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 274-373 of human c-Maf (075444).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075444
Protein Weight	38kDa
Shipping	Ice bag