

C8B Rabbit Polyclonal Antibody For Wb, Elisa - P07358

Item Number: CM0025553



Introduction

C8B Rabbit Polyclonal Antibody targeting complement C8 beta chain, validated for WB and ELISA with cross-reactivity in human, mouse, rat. Ideal for membrane attack complex and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	C8B Rabbit pAb
Remarks/Alias	C82; C8B
Species	Human
GeneID (Human)	732
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-350 of human C8B (NP_000057.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07358
Protein Weight	67kDa
Shipping	Ice bag