

Cadps Rabbit Polyclonal Antibody For Ihc-P, Elisa - Q9Ulu8

Item Number: CM0022499



Introduction

CADPS Rabbit Polyclonal Antibody validated for IHC-P and ELISA. Targets calcium-dependent secretion activator 1 (CAPS-1/CADPS1) in human neural and endocrine tissues. Ideal for exocytosis and neurotransmitter secretion studies.

[Learn More](#)

Parameter	Value
Product Name	CADPS Rabbit pAb
Remarks/Alias	CAPS; CAPS1; CADPS1; UNC-31; CADPS
Species	Human
GeneID (Human)	8618
GeneID	8618
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1020-1100 of human CADPS (NP_899630.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9ULU8
Protein Weight	153kDa
Shipping	Ice bag