

Calpain Small Subunit 1 Rabbit Pab - P04632

Item Number: CM0030886



Introduction

Rabbit polyclonal antibody against Calpain small subunit 1 (CAPNS1). Suitable for Western blot, IHC-P, and ELISA. Detects human, mouse, and rat CAPNS1 protein. Ideal for cytoskeletal remodeling and calpain activity studies.

[Learn More](#)

Parameter	Value
Product Name	Calpain small subunit 1 Rabbit pAb
Remarks/Alias	CANP; CDPS; CSS1; CANPS; CAPN4; CALPAIN4; Calpain small subunit 1
Species	Human
Gene ID (Human)	826
Gene ID	826
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 69-268 of human Calpain small subunit 1 (NP_001740.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04632
Protein Weight	28kDa
Shipping	Ice bag