

Calretinin Rabbit Monoclonal Antibody For Wb, Ihc-P, If-P, Elisa - P22676

Item Number: CM0007602

Introduction

Calretinin Rabbit Monoclonal Antibody (P22676) validated for WB, IHC-P, IF-P, ELISA in Human, Mouse, Rat. Ideal for calcium homeostasis, neuronal plasticity, and brain tissue research.

[Learn More](#)

Parameter	Value
Product Name	Calretinin Rabbit mAb
Remarks/Alias	CR; CAL2; CAB29; Calretinin
Species	Human
Gene ID (Human)	794
Gene ID	794
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-271 of human Calretinin (NP_001731.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22676
Protein Weight	32kDa
Shipping	Ice bag