

Calu Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O43852

Item Number: CM0030100

Introduction

CALU Rabbit Polyclonal Antibody (O43852) validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Calumenin regulates vitamin K-dependent carboxylation and calcium binding. Ideal for ER/Golgi research.

[Learn More](#)

Parameter	Description
Product Name	CALU Rabbit pAb
Remarks/Alias	CALU; calumenin
Species	Human
Gene ID	813
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-315 of human CALU (NP_001210.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43852
Protein Weight	37kDa
Shipping	Ice bag