

Cas9 Rabbit Polyclonal Antibody For Wb And Elisa - Cas9

Item Number: CM0013719

Introduction

CAS9 Rabbit Polyclonal Antibody (pAb) for Western blotting and ELISA applications. Detects CAS9 from Staphylococcus aureus with cross-reactivity to human and other species predicted. Recombinant immunogen. Suitable for CRISPR research.

[Learn More](#)

Parameter	Specification
Product Name	CAS9 Rabbit pAb
Remarks/Alias	CAS9
Species	Staphylococcus aureus
Gene ID (Human)	CAS9
Gene ID	901176
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-730 of Streptococcus pyogenes M1 GAS CAS9 (NP_269215.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Other (Wide Range Predicted)
Shipping	Ice bag