

Casc5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Ng31

Item Number: CM0017026



Introduction

Rabbit polyclonal antibody targeting human CASC5/KNL1, validated for Western blot, IHC-P, and ELISA. Detects human and rat, suitable for kinetochore and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	CASC5 Rabbit pAb
Remarks/Alias	D40; CT29; Spc7; CASC5; MCPH4; hKNL-1; AF15Q14; PPP1R55; hSpc105
Species	Human
Gene ID (Human)	57082
Gene ID	57082
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CASC5 (NP_733468.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8NG31
Protein Weight	265kDa
Shipping	Ice bag