

Casein Kinase 2 Beta (Csnk2B) Rabbit Pab - P67870

Item Number: CM0017743



Introduction

Polyclonal antibody against human CSNK2B, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for CK2 signaling and EBV infection studies.

[Learn More](#)

Parameter	Value
Product Name	Casein Kinase 2 beta (CSNK2B) Rabbit pAb
Remarks/Alias	G5A; CK2B; CK2N; Ckb1; Ckb2; CSK2B; POBINDS; Casein Kinase 2 beta (CSNK2B)
Species	Human
Gene ID (Human)	1460
Gene ID	1460
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-215 of human Casein Kinase 2 beta (Casein Kinase 2 beta (CSNK2B)) (NP_001311.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P67870
Protein Weight	25kDa
Shipping	Ice bag