

Caspase-10 Rabbit Polyclonal Antibody For Wb, Elisa - Q92851

Item Number: CM0013231

Introduction

Caspase-10 Rabbit pAb is a human-specific polyclonal antibody validated for WB and ELISA. It targets the 59 kDa apoptosis initiator caspase, ideal for apoptosis and cell death research.

[Learn More](#)

Parameter	Value
Product Name	Caspase-10 Rabbit pAb
Remarks/Alias	MCH4; ALPS2; FLICE2; FLICE-2; Caspase-10
Species	Human
GeneID (human)	843
GeneID	843
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 220-521 of human Caspase-10 (NP_116759.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q92851
Protein Weight	59kDa
Shipping	Ice bag