

Caspase-12 Rabbit Monoclonal Antibody For Wb, Ihc, If, Elisa - Q6Uxs9

Item Number: CM0007229

Introduction

Rabbit monoclonal antibody targeting human Caspase-12 protein. Suitable for Western blot, IHC, IF, and ELISA applications. Exhibits cross-reactivity with mouse and rat. Useful for innate immunity and inflammatory response research.

[Learn More](#)

Parameter	Value
Product Name	Caspase-12 Rabbit mAb
Remarks/Alias	CASP-12; CASP12P1; Caspase-12
Species	Human
GeneID (Human)	100506742
GeneID	100506742
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human Caspase-12 (NP_001177945.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6UXS9
Protein Weight	39kDa
Shipping	Ice bag