

Caspase-8 Rabbit Pab - O89110

Item Number: CM0013917

Introduction

Rabbit polyclonal antibody to Caspase-8, a key apoptosis, necroptosis, and pyroptosis regulator. Validated for WB, IF/ICC, ELISA; reacts with human and mouse. Ideal for cell death and inflammation research.

[Learn More](#)

Parameter	Specification
Product Name	Caspase-8 Rabbit pAb
Remarks/Alias	MACH; Mch5; FLICE; CASP-8; Caspase-8
Species	Mouse
GeneID (Human)	841
GeneID	12370
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 4-480 of mouse Caspase 8 (NP_033942.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O89110
Protein Weight	55kDa
Shipping	Ice bag