

Cast Rabbit Pab For Wb, Elisa - P20810

Item Number: CM0030989

Introduction

Rabbit polyclonal antibody against human CAST (Calpastatin) for WB and ELISA. Cross-reacts with mouse and rat. Suitable for calpain inhibition, muscle degradation, and protease research.

[Learn More](#)

Parameter	Value
Product Name	CAST Rabbit pAb
Remarks/Alias	BS-17; PLACK; ST
Species	Human
GeneID (Human)	831
GeneID	831
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20810
Protein Weight	77kDa
Shipping	Ice bag