

Castor1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8Wtx7

Item Number: CM0019280



Introduction

CASTOR1 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human CASTOR1, cross-reacts with mouse and rat. Cytosolic arginine sensor for mTORC1 signaling research.

[Learn More](#)

Parameter	Value
Product Name	CASTOR1 Rabbit pAb
Remarks/Alias	GATSL3; CASTOR1
Species	Human
GeneID (Human)	652968
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-329 of human CASTOR1 (NP_001032755.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WTX7
Protein Weight	36kDa
Shipping	Ice bag