

Cathepsin G Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P08311

Item Number: CM0028599



Introduction

Cathepsin G Rabbit pAb validated for WB, IHC-P, IF/ICC, ELISA; reacts with human, mouse, rat. Ideal for studying serine protease functions in inflammation, immune response, and antimicrobial activity.

[Learn More](#)

Parameter	Value
Product Name	Cathepsin G Rabbit pAb
Remarks/Alias	CG; CATG; Cathepsin G
Species	Human
GeneID	1511
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-255 of human Cathepsin G (NP_001902.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08311
Protein Weight	29kDa
Shipping	Ice bag